**EN CONTROL SERUM FOR BIOCHEMISTRY**

Assayed Values for ERBA Reagents in ABNORMAL range  
SysPack Reagents / Liquid Stable Reagents

Cat. No.:	Pack name:	Packaging (Content):
BLT00081	ERBA PATH	R1: 4x5 ml, R2: 1x20 ml

**INTENDED USE**

Control Serum (assayed) is intended for monitoring the accuracy and precision of commonly used clinical chemistry assays in pathological range.

**PRINCIPLE OF THE PROCEDURE**

This human serum control is to be used in the same way as an unknown specimen that would be used in accordance with the direction of the assay procedure. The results obtained for the control are to be compared with the assigned values given in the "assay data" section of the insert and an evaluation made by standard statistical techniques to determine if the procedure is within the control limits.

**COMPOSITION**

The Control serum is prepared from human serum with chemical additives and tissue extracts of human and animal origin. Bacteriostatic agents have been added. The control is provided in lyophilized form for increase stability.

**WARNINGS**

**It is recommended for this product to be handled with same precautions used for patient specimen.**

The blood donations used for production were tested by CE-marked test kits and found to be non-reactive for HBsAg, anti-HIV 1/2 and anti-HCV. In addition HCV and HIV have been tested by PCR. Despite of that the danger of infection for biological material can not be excluded with certainty.

During the work all principles of personal hygiene are to be maintained. Do not eat, drink, smoke. Use the personal protection aids. After finishing the work wash your hands with soap and warm water.

When the bottle is broken the contaminated place must be disinfected.

**RECONSTITUTION**

1. Allow the vial and AQUA-4 (supplied in the kit) to attain room temperature. Add exactly 5 ml of AQUA-4 and allow to stand for 30 minutes at rest in a light-protected place. Swirl occasionally.

2. Swirl the contents gently to ensure homogeneity before using as sample for testing. Avoid formation of foam.

**STORAGE & STABILITY**
**Prior to reconstitution**

The control should be stored at 2–8°C and is stable till the expiry date printed on the label. Protect from light.

**After reconstitution**

The constituents when protected from light and contamination are stable:

- 1 day at 18–25°C
- 1 week at 2–8°C (Bilirubin: 1 day)
- 1 month at –20 °C

**LIMITATIONS**

The results obtained using the control are dependent upon several factors. Erroneous results can occur from reconstitution inaccuracy and the technique errors associated with the assay procedure. The serum is not compatible for use with o-Tolidine Glucose procedures.

Improper storage or handling of the control can also affect the results. If there is a visible evidence of microbial growth in a vial, do not use that vial.

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