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REFERENCES / ЛИТЕРАТУРА / REFERENCIAS

1. Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
2. Barr, D.P., Russ E. M., Eder, H.A., Protein-lipid relationships in human plasma, Am. J. Med., 11;480 (1951).
3. Gordon, T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med., 62;707 (1977).
4. Castelli, W.P. et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55;767 (1977).
5. National Institutes of Health publication No. 93-3095, September, (1993).
6. Williams, P., et al., High density lipoprotein and coronary risk factor, Lancet, 1;72, (1979).
7. Kannel, W.B., Castelli, W.P., Gordon, T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med., 90:85, (1979).
8. Castelli, W. P., et al, Cholesterol and other lipids in coronary heart disease.
9. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 1995; 119:1127
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.

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