ASO CONTROL

Lot: 0000000

Cat. No.	Pack Name	Packaging (Content)
BLT20004	ASO CON	1 x 1 ml

INTENDED USE

Accuracy control for the determination of ASO in serum by turbidimetry and nephelometry.

COMPOSITION

A dilution of human sera containing high levels of ASO with phosphate buffered saline.

Ready to use.

STORAGE AND STABILITY

The expiry date of the product at 2–8 °C is listed on the label.

After first opening the container, the serum can be used for 6 weeks if stored tightly closed at 2–8 °C after use.

Do not freeze.

PRECAUTIONS AND WARNING

- For *in Vitro* Diagnostic use.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.

ASSIGNED VALUES

ASO:

Value	Interval	1 SD	Unit
Titrated against WHO standardisation:			
xxx	xx–xxx	x.xx	IU/ml
Titrated against Siemens standardisation:			
xxx	xx–xxx	x.xx	IU/ml

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