ASO CALIBRATOR

Lot: XXXXXXX

| Cat. No. | Pack Name | Packaging (Content) |
|----------|------------|---------------------|
| XSYS0051 | ASO CAL SH | 1 x 1 ml |

INTENDED USE

Preparation of reference curves for quantitative immunochemical determination of Antistreptolysin (O) (ASO) in human serum.

COMPOSITION

The reagent is prepared by diluting human serum which contains high levels of ASO with phosphate buffered saline containing 1 % bovine serum albumin. Ready to use.

STORAGE AND STABILITY

The expiry date of the product at 2-8 °C is listed on the label. After first opening the container, the serum can be used for 6 weeks if stored tightly closed at 2-8 °C after use. Do not freeze.

PRECAUTIONS AND WARNING

- For *in Vitro* Diagnostic use.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

ASSIGNED VALUES

ASO: xxx IU/ml
xxx IU/ml

WHO standardisation
Siemens standardisation

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