## XL Calibrant A

Cat. No.	Pack Name	Packaging (Content)
XSYS0064	XL CAL A	3 x 500 ml

### INTENDED USE:

Calibrant A is a buffered saline reagent used for one point calibration of Sodium, Potassium and Chloride electrolytes in the ISE module of the XL Clinical Chemistry Auto-analyzer.

**WARNING:** For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit are not classified as dangerous.

**REAGENT PREPARATION:** The reagent is ready to use.

REAGENT COMPOSITION: Active Ingredients	Concentration
Sodium Chloride	140 mmol/l
Potassium Chloride	4 mmol/l
Triton - X 100	< 0.01%
Buffers, Preservatives	< 0.04%

**REAGENT STORAGE & STABILITY:** Store reagents at room temperature (15...25°C). The reagents are stable till the expiry date mentioned on the label.

**PROCEDURE:** The required amount of XL CAL A solution is to be poured into the XL CAL A Bottle of the XL analyzer. Follow instructions for calibration and sample processing as listed in the XL operator's manual.

**PRECAUTIONS:** 1. Avoid contaminating the reagent.  
2. Keep the container tightly closed and store the reagent in the original container.

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