AST/GOT

Cat. No.	Pack Name	Packaging (Content)
XSYS0016	AST/GOT 330	R1: 6 x 44 ml, R2: 3 x 22 ml
XSYS0073	AST/GOT 564 XL-1000	R1: 6 x 72 ml, R2: 6 x 22 ml

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of AST/GOT (Aspartate Aminotransferase) in human serum and plasma.

CLINICAL SIGNIFICANCE

AST/GOT is widely distributed with high concentrations in the heart, liver, skeletal muscle, kidney and erythrocytes. Damage or disease to any of these tissues such as myocardial infarction, viral hepatitis, liver necrosis, cirrhosis and muscular dystrophy may result in raised levels of AST/GOT.

PRINCIPLE

This reagent is based on IFCC recommendations, without pyridoxal phosphate. The series of reactions involved in the assay system is as follows:

1. SGOT / ASAT present in the sample catalyses the transfer of the amino group from L-aspartate to 2-oxoglutarate forming oxaloacetate and L-glutamate.

2. Oxaloacetate in the presence of NADH and Malate dehydrogenase (MDH) is reduced to L-malate. In this reaction NADH is oxidized to NAD⁺. The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due to the oxidation of NADH to NAD⁺.

3. Addition of Lactate dehydrogenase (LDH) to the reagent is necessary to achieve rapid and complete reduction of endogenous pyruvate so that it does not interfere with the assay.

REAGENT COMPOSITION

R1

Tris buffer (pH 7.8)	110 mmol/l
L-aspartic acid	340 mmol/l
LDH	≥ 4000 U/l
MDH	≥ 750 U/l

R2

CAPSO	20 mmol/l
2-oxoglutarate	85 mmol/l
NADH	1.05 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

On board stability: min. 30 days if refrigerated (2–10°C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Loss activity:

within 3 days	at 2–8 °C	< 8 %
within 3 days	at 15–25 °C	< 10 %

Stability

at least 3 months at -20 °C
Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to IFCC formulation (without PDP).

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES ⁴

At 37°C

Men up to 35 U/l

Women up to 31 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	3.84 U/l
Linearity:	390 U/l
Measuring range:	3.84 – 390 U/l

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	103.2	0.60	0.54
Sample 2	313.2	1.68	0.54

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	43.8	0.60	1.37
Sample 2	115.2	1.08	0.92

COMPARISON

A comparison between XL-Systems AST/GOT (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.941 x - 3.96 U/l

r = 0.996

INTERFERENCES

Following substances do not interfere:
bilirubin up to 30 mg/dl, triglycerides up to 2000 mg/dl, haemolysis interferes due to ASAT activity from erythrocytes.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 of the kit contains less than 1.1 % sodium hydroxide which is classified as corrosive substance.

Warning

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent 2 contains less than 0.1% sodium azide which is classified as very toxic and dangerous substance for the environment.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
3. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
4. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.

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