ALKALICKÁ FOSFATASA AMP LIQUID 500 (ALP-AMP L 500)

Kat. č. 10010216

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Cat. No. 10010216

Store (+2 at +8) °C

Reagents for determination of catalytic concentration of alkaline phosphatase in serum and plasma, optimised IFCC method.

Principle

Under alkaline conditions, alkaline phosphatase (alkaline phosphohydrolase of ortho-phosphoric monoesters, EC 3.1.3.1. - ALP), splits 4-nitrophenylphosphate into 4-nitrophenol and phosphate.

The measure of the catalytic concentration of the enzyme is the amount of released 4-nitrophenol, which is evaluated photometrically in a kinetic way.

Reagents

R1 Buffer

2-Amino-2-Methyl-1-Propanol, pH 10.4 (30 °C)	434 mmol/l
Magnesium acetate	2.48 mmol/l
Zinc sulphate	1.24 mmol/l
HEDTA	2.48 mmol/l

R2 Substrate

4-Nitrophenylphosphate	81.6 mmol/l
------------------------	-------------

Composition of mixture

2-Amino-2-Methyl-1-Propanol, pH 10.4 (30 °C)	350 mmol/l
Magnesium acetate	2 mmol/l
Zinc sulphate	1 mmol/l
HEDTA	2 mmol/l
4-Nitrophenylphosphate	16 mmol/l

Preparation and stability of working reagents

Two reagents method – substrate start

Reagents R1 and R2 are liquid, ready to use.

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

After opening, reagents are stable until expiry date at (+2 to +8) °C if stored at appropriate conditions, closed carefully and without any contamination.

Maximum allowable absorbance of the working reagent measured at 420 nm against distilled water is 1.0

Monoreagent method – sample start

Mix 4 portion of reagent R1 with 1 portion of reagent R2.

Stability: 1 week at (+15 to +25) °C in dark
4 weeks at (+2 to +8) °C in dark

Samples

Non – haemolytic serum, plasma (heparin).

Stability of alkaline phosphatase in serum and plasma is 1 week at (+2 to +8) °C.

Calibration

For the calibration, it is recommended to use Lyonorm Calibrator, Cat. No. 10003200.

Quality control

For quality control, it is recommended to use:

Lyonorm HUM N, Cat. No. 10003204

Lyonorm HUM P, Cat. No. 10003206

Procedure

Wavelength	420 (405 – 430) nm
Cuvette	1 cm
Temperature	+37 °C

Serum/reaction mixture ratio 1/51

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 ml	0.800 ml	0.800 ml
Sample	–	–	0.020 ml
Calibrator	–	0.020 ml	–
Distilled water	0.020 ml	–	–
Mix and after 5 min. incubation (at 37 °C) add:			
Reagent 2	0.200 ml	0.200 ml	0.200 ml

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 ml	1.000 ml	1.000 ml
Sample	–	–	0.020 ml
Calibrator	–	0.020 ml	–
Distilled water	0.020 ml	–	–

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Applications for automatic analysers will be supplied on request.

Calculation

$$\text{Alkaline phosphatase } (\mu\text{kat/l}) = \frac{\Delta A_{\text{sample}}/\text{min}}{\Delta A_{\text{cal}}/\text{min}} \times C_{\text{cal}}$$

C_{cal} = calibrator concentration

Unit conversion

1 kat/l = 1 mol/l · s

1 μkat/l = 60 U/l

Reference values

fS ALP (μkat/l) 37 °C

Men

0.67 – 2.15

Women

0.58 – 1.74

Children (males up to 13 years)

< 6.15

Children (females up to 13 years)

< 3.11

Newborns

< 4.17

The range of reference values is only approximate, it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

Performance data

Linearity: up to 17 μkat/l

Limit of detection: 0.06 μkat/l

Low limit of quantification: 0.17 μkat/l

Working range: 0.17 – 17 μkat/l

Precision (at 37 °C)

INTRA-ASSAY n = 20	Mean (μkat/l)	SD (μkat/l)	CV (%)
Sample 1	1.40	0.01	1.00
Sample 2	3.63	0.02	0.44
Sample 3	10.45	0.06	0.58
INTER-ASSAY n = 20	Mean (μkat/l)	SD (μkat/l)	CV (%)
Sample 1	1.41	0.02	1.07
Sample 2	3.62	0.03	0.77
Sample 3	10.46	0.10	0.98

Accuracy

Bias: – 1.3% (1.53 μkat/l)

+ 1.9% (3.74 μkat/l)

Correlation: with commercial available method. Linear regression

N = 100, r = 0.999, y = 0.928x – 0.083 μkat/l

Interference

Following substances do not interfere:

bilirubin up to 1000 μmol/l, triglycerides up to 20 mmol/l, hemoglobin up to 7 g/l

Note

- If catalytic concentration of alkaline phosphatase exceeds 17 μkat/l (ΔA/min > 0.240), dilute the sample in 1+9 ration with saline (0.9 % NaCl). Multiply the results by 10.
- Procedure for calibration with use of molar absorbance coefficient is available on request.

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagent 1 is not classified but it contains sodium azide (<0.1%) - very toxic substance, hazardous to environment.

Reagent 2 contain sodium azide (<0.1%) - very toxic substance, hazardous to environment, potassium hydroxide (<0.7%) - corrosive substance and acid hydrochloric (<0.4%).

Warning

Hazard statements:

H315 Causes skin irritation

H319 Causes serious eye irritation

Precautionary statements:

P280 Wear protective gloves/protective clothing/eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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