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(AST UV 10x100 P)

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Cat. No. 10003290

Store at (+2 to +8)°C

A set of reagents for preparation of 10x100 ml (AST UV 10x100 P) of working solution of determination of the catalytic concentration of AST in serum and plasma in UV area. Method recommended by DGKCh "new" 1994 (IFCC at 37°C).

Principle

L-aspartate + 2-oxoglutarate $\xrightarrow{\text{AST}}$ oxalacetate + L-glutamate
oxalacetate + NADH + H⁺ $\xrightarrow{\text{MD}}$ L-malate + NAD⁺

AST catalytic concentration is proportional to absorbance drop at 340, 334 or 365 nm.

Linearity

The method is linear up to 7.0 µkat/l for the serum/reaction mixture volume ratio 1/12 (manual assay).
The method is linear up to 10 µkat/l for the serum/reaction mixture volume ratio 1/14.75 (analyser).

References

International Federation of Clinical Chemistry, Scientific Committee: Clin. Chem. 18, 521 (1980)

Reagents

- Enzyme – NADH (10 vials)
LD ≥ 1.0 µkat, MD ≥ 0.7 µkat,
NADH ≥ 21.6 µmol/vial
- Starter (4x28 ml)
2-oxoglutarate 166 mmol/l
Sodium azide 0.1 %
- Buffer solution–substrate (5x200 ml)
Tris 100 mmol/l,
L-aspartate 350 mmol/l
Sodium azide 0.1 %
- Activator (20 tbl)
Pyridoxal-5-phosphate 6 µmol/tablet

Reaction mixture

Tris buffer, pH 7.8 (37°C) 80 mmol/l
L-aspartate 290 mmol/l
Lactatdehydrogenase (LD) ≥ 10 µkat/l
Malatdehydrogenase (MD) ≥ 7 µkat/l
NADH ≥ 180 µmol/l
2-oxoglutarate 13.8 mmol/l
pyridoxal-5-fosphate 100 µmol/l
Serum/incubation mixture volume ratio 1/12

Reference values

fS AST (µkat/l) 37°C
men 0.17–0.85
women 0.17–0.60
The range of reference values is only approximate, it is recommended to all laboratories to verify the extension of reference interval for their concrete examined population.

Calibration and quality control

Reference attested sera with stated catalytic concentration for AST using the determined stated method. We recommend: BIO-LA-TEST® LYONORM CALIBRATOR, Cat. No. 10003200
BIO-LA-TEST® LYONORM HUM N, Cat. No. 10003204
BIO-LA-TEST® LYONORM HUM P, Cat. No. 10003206

Preparation of working solution

Dissolve content of the vial with the Reagent 1 in 100 ml of Reagent 3. After dissolution add 2 tablets of Reagent 4.

Stability: 14 days at (+2 to +8)°C in a dark
2 days at (+15 to +25)°C in a dark

Procedure

Samples: non-haemolytic serum, heparinized or EDTA plasma
Wavelength 340, 334, 365 nm
Cuvette 1 cm
Temperature (37 ± 0.1)°C

Measure (ml)	Sample
Working solution	1.00
Serum	0.10
Stir, incubate at 37 °C for 10 minutes and add	
Reagent 2	0.10
Stir, incubate at 37 °C for 2 minutes and measure absorbance in 1- min. intervals for at least 3 minutes. Calculate absorbance average deviation within 1 min. (ΔA).	

Calculation

ALT (µkat/l) = ΔA x 31.7 (340 nm)
ΔA x 32.4 (334 nm)
ΔA x 58.8 (365 nm)

Notes

We don't recommend to freeze the samples before the analysis, they should be only stored at (+2 to +8) °C for 4 days as max.

If during manual assay sample value exceeds 7.0 µkat/l, and during assay on an analyser it is up to 10 µkat/l, the sample must be diluted with physiological saline (NaCl) in the proportion 1 + 4 (multiply result 5x).

The analysis may be carried out at 30°C without substantial change in the method characteristics. The following values are modified:

Reference values:

fS AST (µkat/l) 30°C men 0.1–0.5
women 0.1–0.2

Conversion of units

1 kat/l = 1 mol/l.s
1 µkat/l = 60 U/l

Stability and storage

Store reagent at (+2 to +8)°C before and after opening in dark.

Expiration time is mentioned on the cover.

After using the vial close immediately to prevent the evaporation or contamination of the reagent.

Performance characteristics

Linearity: up to 10 µkat/l
Detection limit: 0.039 µkat/l
Measuring range: 0.12 - 10 µkat/l

Precision (at 37°C)

INTRA-ASSAY n =20	Mean (µkat/l)	SD (µkat/l)	CV (%)
Sample 1	0.84	0.04	4.60
Sample 2	2.69	0.03	1.25
Sample 3	4.39	0.06	1.29
INTER-ASSAY n =10	Mean (µkat/l)	SD (µkat/l)	CV (%)
Sample 1	0.85	0.04	4.76
Sample 2	2.71	0.04	1.35
Sample 3	4.22	0.07	1.61

Accuracy

Test was compared with commercially available test. Results:

N = 100, r = 0.998, y = 1.028 x -0.0405 µkat/l

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

The set contains sodium azide in low concentration (0.1 %) in reagents 2 and 3, classified as very toxic and dangerous substance to the environment.

Reagent 2

- R 22 Harmful if swallowed.
S 28 After contact with skin, wash immediately with plenty of water.
S 45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Reagent 3

- R 22 Harmful if swallowed.
R 36/38 Irritating to eyes and skin.
S 28 After contact with skin, wash immediately with plenty of water.
S 45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

S 36/37/39 Wear suitable protective clothing, gloves and eye/face protection.

In the case of accidental ingestion rinse your mouth and drink about 0.5 l water, in the case of skin staining or affection of eyes rinse the skin or eyes with a stream of clean water. Seek for medical attention.

Waste disposal

Any processed sample should be treated as a potentially infectious one and together with possible residues of reagents they should be disposed in accordance with the company internal regulations as dangerous waste in compliance with the Waste Act.

Paper and other rinsed packages should be disposed depending on kind of the material as sorted waste (paper, glass, plastic).

Date of last revision: 1. 4. 2011