GLUCOSE

Cat. No.	Pack Name	Packaging (Content)
BLT00025	GLU 500	R1: 2 x 250 ml, R2 standard: 2 x 5 ml
BLT00026	GLU 4x250	R1: 4 x 250 ml, R2 standard: 2 x 5 ml
BLT00027	GLU 1000	R1: 1 x 1000 ml

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Glucose in human serum, plasma and urine.

CLINICAL SIGNIFICANCE



Accurate measurement of glucose in body fluid is important in diagnosis and management of diabetes, hypoglycemia, adrenal dysfunction and various other conditions. High levels of serum glucose may be seen in case of Diabetes mellitus, in patients receiving glucose containing fluids intravenously, during severe stress and in cerebrovascular accidents.

Decreased levels of glucose can be due to insulin administration, as a result of insulinoma, inborn errors of carbohydrate metabolism or fasting.

PRINCIPLE

Trinder's method

Glucose in the sample is oxidised to yield gluconic acid and hydrogen peroxide in the presence of Glucose oxidase. The enzyme peroxidase catalyses the oxidative coupling of 4-aminoantipyrine with phenol to yield a coloured quinonimine complex, with absorbance proportional to the concentration of glucose in sample.

REAGENT COMPOSITION

R1

Phosphate buffer	250 mmol/l
Glucose oxidase	>25 U/ml
Peroxidase	>2 U/ml
Phenol	5 mmol/l
4-aminoantipyrine	0.5 mmol/l

R2 standard	See bottle label
--------------------	------------------

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability after addition of a glycolytic inhibitor (Fluoride, monoiodoacetate, mannose): ⁴

2 days	at 20–25 °C
7 days	at 4–8 °C

Stability in serum (separated from cellular contents, hemolysis free) without adding a glycolytic inhibitor: ^{2,5}

8 hours	at 25 °C
72 hours	at 4 °C

Stability of glucose in urine: 24 hours at 4–8 °C

For the determination in urine dilute the sample using redistilled water in 1 + 10 ratio. Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit or the calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

mg/dl x 0.056 = mmol/l

EXPECTED VALUES ²

Serum:

Glucose Fasting:

Cord:	45 – 96 mg/dl
Newborn, 1 d:	40 – 60 mg/dl
Newborn, >1 d:	50 – 80 mg/dl
Child:	60 – 100 mg/dl
Adult:	74 – 100 mg/dl
> 60 y:	82 – 115 mg/dl
> 90 y:	75 – 121 mg/dl

Glucose 2 h Postprandial:

<120 mg/dl

WB (Hep) Adult: 65 – 95 mg/dl

Urine: 1 – 15 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 2.34 mg/dl

Linearity: 450 mg/dl

Measuring range: 2.34 – 450 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	153.75	1.61	1.05
Sample 2	239.64	3.93	1.66

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	97.14	1.61	1.64
Sample 2	258.39	3.04	1.19

COMPARISON

A comparison between XL-Systems Glucose (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.000 x - 0.714 \text{ mg/dl}$$

$$r = 0.999$$

INTERFERENCES

Following substances do not interfere:
haemoglobin up to 7.5 g/l, bilirubin up to 30 mg/dl, triglycerides up to 750 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent R1 contains 0.74% NaOH and is classified as irritant and also contains less than 0.1% sodium azide which is classified as very toxic and dangerous substance for the environment.

Warning

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength (490 – 550) maximum at 500 nm

Cuvette 1 cm

	Reagent blank	Standard (Cal.)	Sample
Reagent 1	1000 µl	1000 µl	1000 µl
Sample	-	-	10 µl
Standard (Cal.)	-	10 µl	-
Distilled water	10 µl	-	-

Mix and incubate 5 - 10 min. at 37 °C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank.

CALCULATION

$$\text{Glucose (mg/dl)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{st}}} \times C_{\text{st}} \quad C_{\text{st}} = \text{standard (calibrator) concentration}$$

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	505
Wavelength 2 (nm)	670
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	5
Incubation temp. (°C)	37
Normal Low (mg/dl)	74
Normal High (mg/dl)	100
Linearity Low (mg/dl)	2.34
Linearity High (mg/dl)	450
Concentration of Standard	See bottle label
Blank with	Reagent
Absorbance limit (max.)	0.2
Units	mg/dl

Erba Lachema s.r.o., Karásek 1d, 621 00 Brno, CZ
e-mail: diagnostics@erbalachema.com, www.erbamannheim.com

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

REFERENCES / ЛИТЕРАТУРА / LITERATURA / LITERATŪRA / REFERENCIAS

1. Thomas L.: Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 131 - 7.
2. Tietz N. W., (Ed.), Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
3. Barham, D., Trinder, P.: An improved color reagent for the determination of blood glucose by the oxidase system. Analyst, 1972, 97; 142 - 5.
4. Guder WG, Zawta B et al. The quality of Diagnostic Samples. 1st ed. Darmstadt: GIT verlag; 2001;p.30-1.
5. Snacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chemi 2002; 48:436-72.

SYMBOLS USED ON LABELS / СИМВОЛЫ, ИСПОЛЬЗУЕМЫЕ НА ЭТИКЕТКАХ / SYMBOLY, POUŽITÉ NA ETIKETÁCH / SÍMBOLOS UTILIZADOS EN LAS ETIQUETAS

<p>REF Catalogue Number Katalogové číslo Katalogové číslo Номер каталога Code de Catalogue Código de Catalogo</p>	<p> Manufacturer Výrobce Výrobca Производитель Fabriqué par... Fabricado por...</p>	<p> See Instruction for Use Čtěte návod k použití Čítajte návod k použitiu Перед использованием Внимательно изучайте инструкцию Lire les Instructions avant l'Utilisation Ver Instrucciones Para su Uso</p>
<p>LOT Lot Number Číslo šarže Номер партии Número de Lot Número de Lote</p>	<p> CE Mark - Device comply with the Directive 98/79/EC Знак CE - соответствие Директиве 98/79/EC</p>	<p> Storage Temperature Teplota skladování Teplota skladovania Температура хранения Limites de Température Rango de Temperatura</p>
<p> Expiry Date Datum expirace Dátum expirácie Срок годности Date d'Expiration Fecha de Vencimiento</p>	<p>IVD In Vitro Diagnostics In vitro Diagnostikum Ин витро диагностика Trousse Médicale Diagnostique in Vitro Dispositivo Médico para Diagnóstico in Vitro Solamente</p>	<p>CONT Content / Obsah / Содержание Contenus / Contenido</p> <p> Национальный знак соответствия для Украины Ukrainian quality mark</p>

QUALITY SYSTEM CERTIFIED
ISO 9001 ISO 13485

 Erba Lachema s.r.o., Karásek 1d, 621 00 Brno, CZ
e-mail: diagnostics@erbalachema.com, www.erbamannheim.com

N/55/13/D/INT Date of revision: 4. 3. 2015