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(Cu 35)

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(Cu 35)

Cat. No 10003162 Store at (+2 to +8)°C

Reagent Set for the preparation of 100 ml of working solutions used for the determination of copper in serum. The Set is sufficient for at least 35 analyses.

Principle
Cu(II) ions form with Bathocuproine a stable orange coloured complex.

References
Landers, J.W., Zak, B.: J. Clin. Path. 29, 590 (1958)
Szilágyi, L., Páhoki, I.: Orv. Hetil. 109, 2837 (1968)

Reagents

1	Standard solution	(2 ml)
	copper sulfate 3 mmol/l	
2	Bathocuproine	(50 ml)
	solution of 0.44 mmol/l of 2,9-dimethyl-4,7-diphenyl-1,10-phenanthroline-3,6-disulfonate (sodium salt), sodium acetate 2 mol/l	
3	Reducing agent	(1 vial)
	hydroxylamine sulfate 1.25 g/vial	
4	Deproteinizing agent	(50 ml)
	trichloroacetic acid 0.6 mol/l, hydrochloric acid 2 mol/l	

Composition of reaction mixture

Bathocuproine	0.22 mmol/l
Trichloroacetic acid	0.15 mol/l
Hydrochloric acid	0.50 mol/l
Sodium acetate	1.00 mol/l
Hydroxylamine sulfate	0.04 mol/l
Volume ratio of serum/reaction mixture	1/2

Reference values
fS Copper (µmol/l)

Men	11.0 – 22.0
Women	13.4 – 24.4

The range of reference values is only approximate, it is recommended to all laboratories to verify the extension of reference interval for their concrete examined population.

Quality control
For control it is recommended to use control materials with values for values for copper, for example Bio Rad Lyphochek Assayed Chemistry control.

Working solutions

- Solution 1** Into 100 ml volumetric flask pipette 1 ml of Reagent 1 and dilute with distilled water up to the mark. Solution contains 30 µmol of Cu(II)/l.
Stability: several months when stored at (+15 to +25)°C.
- Solution 2** Dissolve the content of the vial with Reagent 3 in the entire volume of Reagent 4.
Stability: 3 months when stored at (+2 to +8)°C.

Procedure

Samples: serum
Wavelength: (460 – 500) nm
Cuvette: 1 cm
Temperature: (+15 to +25)°C

In three test tubes mix Solution 2 in the ratio 1+1 with serum (sample), Solution 1 (standard), or with distilled water (control solution). Stir and after 30 min centrifuge the sample for 10 min at 3000 rpm (for example, mix 1 ml of Solution 2 and 1 ml of serum or Solution 1, or distilled water). In another three test tubes mix supernatant, standard, or control solution in the ratio 1+1 with Reagent 2 (for example, mix 1 ml of supernatant, standard, or Control solution with 1 ml of Reagent 2). Stir and between (5 to 90) min read the absorbances of sample (A_1) and standard (A_2) against control solution.

Pipette (ml)	Sample	Standard	Control soln.
Solution 2	1.00	1.00	1.00
Serum	1.00	–	–
Solution 1	–	1.00	–
Dist. water	–	–	1.00
Mix, let stand for 30 min, centrifuge for 10 min at 3000 rpm and pipette into another test tube			
Supernatant	1.00	1.00	1.00
Reagent 2	1.00	1.00	1.00
Mix and within (5–90) min read the absorbances of the sample (A_1) and standard (A_2) against control solution.			

Calculation

$$\text{Copper } (\mu\text{mol/l}) = 30 \cdot \frac{A_1}{A_2}$$

Notes

The laboratory glassware for the determination of copper must be perfectly clean and used for this purpose only. After washing the glassware by the ordinary way, it is recommended to dip it overnight into approx. 2 % solution of EDTA in ammonia diluted in the ratio 1+1 with distilled water, then rinse it thoroughly with redistilled water and dry well.

Performance characteristic

Linearity: up to 100 µmol/l
Detection limit: 1.04 µmol/l
Measuring range: 3.46 – 100 µmol/l

Precision (at 37°C)

INTRA-ASSAY	Mean (µmol/l)	SD (µmol/l)	CV (%)
n = 20			
Sample 1	18.8	0.73	3.86
Sample 2	35.4	1.42	4.00
Sample 3	42.6	1.20	2.82
INTER-ASSAY	Mean (µmol/l)	SD (µmol/l)	CV (%)
n = 10			
Sample 1	19.1	0.85	4.46
Sample 2	36.1	1.25	4.34
Sample 3	43.3	1.56	3.45

Accuracy

Test was compared with commercially available test. Results:
N = 26, r = 0.965, y = 1.245 x - 2.582 µmol/l

Health protection

For *in vitro* diagnostic use.
To be handled by entitled and professionally educated person.

Reagent 3 contains hydroxylamine sulfate (100 %) which is toxic (category 4), irritant, carcinogenic and toxic substance for environment.

Warning

Hazard statement:

H290 May be corrosive to metals.
H302 Harmful if swallowed.
H312 Harmful in contact with skin.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H319 Causes serious eye irritation.
H351 Suspected of causing cancer.
H373 May cause damage to organs through prolonged or repeated exposure.
H400 Very toxic to aquatic life.

Precautionary statement:

P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection.
P302+P352 IF ON SKIN: Wash with plenty of water and soap.
P308+P313 IF exposed or concerned: Get medical advice.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent 4 contains trichloroacetic acid (9 %) - corrosive substance and toxic substance for environment and hydrochloric acid (18 %) - corrosive substance.

Danger

Hazard statement:

H314 Causes severe skin burns and eye damage.
H335 May cause respiratory irritation.
H411 Toxic to aquatic life with long lasting effects.

Precautionary statement:

P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of water and soap.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

First aid

At an accidental ingestion drink of approx. 0.5 ml of water. At an eye contact flush immediately and thoroughly with large quantity of water. At a skin contact wash the affected spot with soap and warm water. In all serious cases of health damage consult the physician.

Stability and storage

Store the set at (+2 to +8)°C. The expiration date is stated on the package.
Close the vial immediately after using so that the evaporation or contamination of the reagent does not occur. Store at recommended temperature and in the dark.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging paper waste and rinsed containers to recycling.

Date of last revision: 23. 2. 2015