RU

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CE

IVD

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(Fe 70)

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  BIO-LA-TEST®  																	
IRON (Fe 70)																	
Cat. No. 10003155	Store at (+2 to +8)°C																
Set of reagents for preparation of 210 ml of working solutions for iron determination in serum. The reagents are sufficient for 70 analyses at least. The set can be used also for determination of total binding capacity in serum, see working instruction of BLT Total binding capacity.																	
Principle Bathophenanthroline forms a red complex with ferrous ions. This complex is suitable for photometric determination.																	
References International Committee for Standardization in Haematology: Clin. Chem. 17, 127 (1971) Caraway, W. T.: Clin. Chem. 9, 188 (1963) Zak, D.: Clin. Chim. Acta 24, 293 (1958)																	
Reagents <table border="0" style="width: 100%;"> <tr> <td>1 Standard solution</td> <td style="text-align: right;">(40 ml)</td> </tr> <tr> <td colspan="2">ferrous ammonium sulfate 17.9 µmol/l</td> </tr> <tr> <td>2 Bathophenanthroline</td> <td style="text-align: right;">(110 ml)</td> </tr> <tr> <td colspan="2">4.7-diphenyl-1.10-phenanthroline-3.6-disulphonic acid 0.46 mmol/l,</td> </tr> <tr> <td colspan="2">sodium acetate 2 mol/l</td> </tr> <tr> <td>3 Thioglycolic acid 92.1 mol/l</td> <td style="text-align: right;">(4 ml)</td> </tr> <tr> <td>4 Deproteinizing reagent</td> <td style="text-align: right;">(110 ml)</td> </tr> <tr> <td colspan="2">trichloroacetic acid 0.6 mol/l, hydrochloric acid 2 mol/l</td> </tr> </table>		1 Standard solution	(40 ml)	ferrous ammonium sulfate 17.9 µmol/l		2 Bathophenanthroline	(110 ml)	4.7-diphenyl-1.10-phenanthroline-3.6-disulphonic acid 0.46 mmol/l,		sodium acetate 2 mol/l		3 Thioglycolic acid 92.1 mol/l	(4 ml)	4 Deproteinizing reagent	(110 ml)	trichloroacetic acid 0.6 mol/l, hydrochloric acid 2 mol/l	
1 Standard solution	(40 ml)																
ferrous ammonium sulfate 17.9 µmol/l																	
2 Bathophenanthroline	(110 ml)																
4.7-diphenyl-1.10-phenanthroline-3.6-disulphonic acid 0.46 mmol/l,																	
sodium acetate 2 mol/l																	
3 Thioglycolic acid 92.1 mol/l	(4 ml)																
4 Deproteinizing reagent	(110 ml)																
trichloroacetic acid 0.6 mol/l, hydrochloric acid 2 mol/l																	
Reaction mixture <table border="0" style="width: 100%;"> <tr> <td>Bathophenanthroline</td> <td style="text-align: right;">0.23 mmol/l</td> </tr> <tr> <td>Thioglycolic acid</td> <td style="text-align: right;">95.50 mmol/l</td> </tr> <tr> <td>Trichloroacetic acid</td> <td style="text-align: right;">0.15 mol/l</td> </tr> <tr> <td>Hydrochloric acid</td> <td style="text-align: right;">0.50 mol/l</td> </tr> <tr> <td>Sodium acetate</td> <td style="text-align: right;">1.00 mol/l</td> </tr> <tr> <td>Volume ratio serum/reaction mixture</td> <td style="text-align: right;">1/2</td> </tr> </table>		Bathophenanthroline	0.23 mmol/l	Thioglycolic acid	95.50 mmol/l	Trichloroacetic acid	0.15 mol/l	Hydrochloric acid	0.50 mol/l	Sodium acetate	1.00 mol/l	Volume ratio serum/reaction mixture	1/2				
Bathophenanthroline	0.23 mmol/l																
Thioglycolic acid	95.50 mmol/l																
Trichloroacetic acid	0.15 mol/l																
Hydrochloric acid	0.50 mol/l																
Sodium acetate	1.00 mol/l																
Volume ratio serum/reaction mixture	1/2																
Reference values fS Iron (µmol/l) men 14.3–26.0 women 10.7–21.5 The range of reference values is only approximate, it is recommended to all laboratories to verify the extension of reference interval for their concrete examined population.																	
Quality control BIO-LA-TEST® LYONORM U, Cat. No. 10003176 BIO-LA-TEST® LYONORM P, Cat. No. 10003172 Control human sera with iron determined as mentioned above are recommended.																	

Preparation of working solution

Into a 100 ml graduated flask, pipette 3 ml of Reagent 3. Adjust to 100 ml with Reagent 4. In case of precipitate formation filter through close filter.

Stability: 2 months under +8°C in the dark.

Procedure

Samples:	non-haemolytic serum
Wavelength	(510–550) nm
Cuvette	1 cm
Temperature	(+15 to +25)°C

In one test tube mix together 1.00 ml of working solution and 1.00 ml of either serum or VK¹ (sample). In other two test tubes pipette 0.50 ml of working solution and 0.50 ml of Reagent 1 (standard) and 0.50 ml of working solution and 0.50 ml of bidistilled water, respectively (blank). Mix thoroughly and after 5 minutes at (+15 to +25)°C centrifuge the sample tube 10 min at 3000 rpm.

Into a test tube pipette 1.00 ml of the supernatant. Into this tube as well as to the standard and control tubes add 1.00 ml of Reagent 2, and mix. Between 5 and 60 min read the absorbances of the sample (A₁) and standard (A₂) against the blank.

Calculation

$$\text{Iron } (\mu\text{mol/l}) = 17.90 \cdot \frac{A_1}{A_2}$$

Notes

VK¹ is the supernatant used in total binding capacity determination using BIO-LA-TEST® TOTAL BINDING CAPACITY, Cat. No. 10003154, which is successor of BIO-LA-TEST® IRON.

The glassware used for determination of iron must be perfectly clean, preferably assigned to this method only. After customary washing it is advisable to soak the glassware overnight in about 2% Na₂EDTA in diluted ammonia 1+1, than thoroughly rinse with bidistilled water and dry. Light pink coloration of the reagent 2 is acceptable.

Performance characteristic

Linearity: up to 150 µmol/l

Detection limit: 1.05 µmol/l

Measuring range: 3.51 - 150 µmol/l

Precision (at 37°C)

INTRA-ASSAY n=20	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	21.8	0.57	2.62
Sample 2	30.3	0.97	3.20
Sample 3	33.4	1.08	3.24
INTER-ASSAY n=10	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	20.9	0.87	4.16
Sample 2	29.5	1.00	3.39
Sample 3	32.1	1.22	3.78

Accuracy

test was compared with commercially available test. Results:
 N=44, r=0.983, y=0.82 x -4.39 µmol/l

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Reagent 3 contains thioglycolic acid (100%) which is toxic and corrosive. Reagent 4 contains hydrochloric acid (7.2%) and trichloroacetic acid (10%) which are corrosive.

The working solution contains all these reagents.

All the principles of personal hygiene must be observed. Do not eat, drink or smoke when working. Wash your hands thoroughly after work.

On skin or eye contact rinse the contact place with excess of tap water. When swallowed, drink 0.5 l of water. On vapour intoxication provide fresh air. Consult a physician!

Waste disposal

All tested samples should be treated as potentially infectious and with the containing rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and rinsed containers to recycling.

Date of last revision: 21. 1. 2011

Produced Erba Lachema s.r.o. Karásek 1d, 62133 Brno, CZ e-mail: diagnostics@lachema.com , www.lachema.com
--