Železo Liquid 200

(Fe L 200)

Kat.č. 10007190

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Cat. No. 10007190

Store at (+2 to +8) °C

A set of reagents for direct colorimetric determination of iron without deproteination in serum and plasma.

Principle

In a pH 4.8 buffer system, iron is released from transferrin and then quantitatively reduced to ferrous state. Fe²⁺ forms with Ferene S ((3-(2-pyridil)-5,6-bis-2-(5-furysulfonic acid)-1,2,3-triazine)) a stable coloured complex, which colour intensity is proportional to the amount of iron in the sample. The interference from copper is eliminated by particular reaction conditions and a specific masking agent.

Reagents

R1 Reagent 1

Acetate buffer pH 4.8
Guanidine hydrochloride
Thiourea

4 x 45 ml

1.4 mol/l
≥ 4.5 mol/l
65 mmol/l

R2 Reagent 2

Ferene S
Ascorbic Acid

1 x 18 ml
≥ 20 mmol/l
≥ 500 mmol/l

R3 Standard

Ferrous ammonium sulphate

1 x 40 ml
17.9 mmol/l

Reaction mixture

Acetate buffer pH 4.8
Guanidine hydrochloride
Thiourea
Ferene S
Ascorbic Acid

1.12 mol/l
≥ 3.6 mol/l
52 mmol/l
≥ 8 mmol/l
≥ 20 mmol/l

Preparation and stability of working solutions

If store at (+2 to +8) °C, reagents in unopened vials are stable up to expiry date indicated on package.

Reagents are liquid and ready to use. After the first opening, reagents are stable 90 days at (+2 to +8) °C, if contamination avoided and vials are recapped immediately after use.

Standard is liquid and ready to use.

Samples

Serum or plasma not haemolyzed. Use only heparin salts as anticoagulants. For TIBC determination use a supernatant received by a procedure described at the working instruction for BLT Total binding capacity. Stability of the samples: 7 days at (+2 to +8) °C.

Calibration

For calibration it is recommended:

LYONORM CALIBRATOR, Cat. No. BLT00069 or the standard R3 included in the set.

Quality control

For control it is recommended:

LYONORM HUM N, Cat. No. BLT00070
LYONORM HUM P, Cat. No. BLT00071

Procedure

Wavelength 593 (580–600) nm
Cuvette 1 cm
Temperature 37 °C
Volume ration serum/reaction mixture 1/6,25
Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

	Reagent blank	Sample blank	Sample	Standard (calibrator)
Distilled water	0.2 ml	0.05 ml	–	–
Sample	–	0.2 ml	0.2 ml	–
Standard (calibrator)	–	–	–	0.2 ml
Reagent R1	1.0 ml	1.0 ml	1.0 ml	1.0 ml
Reagent R2	0.05 ml	–	0.05 ml	0.05 ml

Mix and after 5 minutes incubation at 37 °C, read the absorbance of sample blank A1 against distilled water and the absorbance of sample A2 and standard (calibrator) A3 against reagent blank.

Applications for automatic analysers will be supplies on request.

Calculation

$$\text{Iron } (\mu\text{mol/l}) = \frac{A_2 - A_1}{A_3} \times C_{st}$$

C_{st} = standard (calibrator) concentration

Conversion of units

1 μmol/l = 0.179 μg/dl

Reference values

fS iron (μmol/l)

male: 10.7–28.6
female: 7.2–25.9

The range of reference values is only approximate, it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

Performance data

Linearity: up to 240 μmol/l. Samples with higher concentrations dilute 1:10 with normal saline (for example 1 ml of sample + 9 ml of saline) and results multiplied by 10.

Limit of detection: 0.23 μmol/l

Low limit of quantification: 0.7 μmol/l

Working range: 0.7–240 μmol/l

Precision (at 37 °C)

INTRA-ASSAY n = 20	Mean (μmol/l)	SD (μmol/l)	CV (%)
Sample 1	9.76	0.08	0.79
Sample 2	19.09	0.14	0.75
Sample 3	33.89	0.22	0.65
INTER-ASSAY n = 10	Mean (μmol/l)	SD (μmol/l)	CV (%)
Sample 1	8.78	0.13	1.34
Sample 2	19.51	0.20	1.02
Sample 3	33.99	0.29	0.86

Accuracy: (control materials validated in EQA)

Bias: +0.35% (18.4 μmol/l)

–0.32% (36.2 μmol/l)

Correlation: with commercial available method. Linear regression:

N = 100, y = 1.017 x + 0.602 μmol/l, r = 0.999, s_{x,y} = 0.797

Interferences: bilirubin up to 256 μmol/l, haemoglobin up to 0.5 g/l and triglycerides up to 11.3 mmol/l do not interfere.

Note

- It is recommended to use disposable test tubes for iron analysis. In case you use glassware, it must be perfectly clean and used for this purpose only. It is recommended to wash it in a solution of 1mol/l hydrochloric acid and distilled water. After washing it is recommended to dip glassware overnight into about 2% EDTA solution in 10% ammonia, then rinse well with redistilled or deionized water and dry.
- Reagent R2 is a yellow-colored solution that may turn slightly darker, such a change will not effect the reagent performance.
- The both reagents R1 a R2 must be without any turbidity, use only transparent solutions.

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent R1 contains guanidine hydrochloride (38%), acetic acid (7%) and sodium hydroxide (3%).

Danger

Hazard statement:

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P301+P312 IF SWALLOWED: Call a doctor if you feel unwell.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

First aid

In case of an accidental ingestion wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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