# PHOSPHORUS

Cat. No.	Pack Name	Packaging (Content)
BLT00047	PHOS 100	R1: 2 x 50 ml, R2 standard: 1 x 5 ml
BLT00048	PHOS 250	R1: 1 x 250 ml, R2 standard: 1 x 5 ml

EN

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Phosphorus in human serum, plasma or urine.

## CLINICAL SIGNIFICANCE

More than 80% of the body's phosphate is present in bones as calcium phosphate. The remainder is found intracellularly as organic phosphates such as phospholipids, nucleic acids and ATP or extracellularly as inorganic phosphorus.

There is generally a reciprocal relationship between serum calcium and inorganic phosphorus levels. Increased levels of serum phosphorus is seen in renal diseases, hypoparathyroidism and excessive vitamin D intake.

Decreased levels of phosphorus is seen in rickets, osteomalacia (adult rickets), hyperparathyroidism and in diabetic coma.

## PRINCIPLE

Inorganic phosphorus combines with ammonium molybdate in the presence of strong acids to form phosphomolybdate. The formation of reduced phosphomolybdate is measured at 340 nm and is directly proportional to the concentration of inorganic phosphorus present in the sample.

## REAGENT COMPOSITION

### R1

Ammonium Molybdate 0.35 mmol/l  
Sulphuric acid 131 mmol/l

**R2 standard** See bottle label

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

## SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

### Stability in serum / plasma:

7 days at 4–25°C  
3 months at -20°C

### Stability in urine:

2 days at 20–25°C at pH < 5

For collection of 24 h urine add 20 to 30 ml of 6 mol/l HCl into the collection bottle to avoid phosphate precipitations.

Dilute urine 1 + 10 with dist. water before determination and multiply the result by 11. Discard contaminated specimens.

## CALIBRATION

Calibration with the standard included in the kit or the calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

## QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

## UNIT CONVERSION

mg/dl x 0.32 = mmol/l

## EXPECTED VALUES †

### Serum

Adult 2.5 – 4.5 mg/dl  
Children 4.0 – 7.0 mg/dl

### Urine, 24 h

Adult 0.4 – 1.3 g / 24 h

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.17 mg/dl

**Linearity:** 23.25 mg/dl

**Measuring range:** 0.17 – 23.25 mg/dl

## PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	4.15	0.10	2.43
Sample 2	3.90	0.05	1.30

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	3.99	0.15	3.75
Sample 2	11.63	0.29	2.51

## COMPARISON

A comparison between XL-Systems Phosphorus (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.982 x - 0.118 mg/dl

r = 0.992

## INTERFERENCES

Following substances do not interfere:  
haemoglobin up to 1.25 g/l, bilirubin up to 20 mg/dl, triglycerides up to 500 mg/dl.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous but R2 standard contains less than 0.1% sodium azide - classified as very toxic and dangerous substance for the environment.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## ASSAY PROCEDURE

**Wavelength:** 340 nm, Hg 334 nm, Hg 365 nm

**Cuvette:** 1 cm

	Reagent blank	Standard (Cal.)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	-	-	0.01 ml
Standard (Cal.)	-	0.01 ml	-
Distilled water	0.01 ml	-	-

Mix and incubate 5 min. at 37°C. Measure absorbance of the sample  $A_{sam}$  and standard  $A_{st}$  against reagent blank.

## CALCULATION

Phosphorus (mmol/l) =  $\frac{\Delta A_{sam}}{\Delta A_{st}} \times C_{st}$       $C_{st}$  = standard (calibrator) concentration

Applications for automatic analysers are available on request.

## ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	340
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	5
Incubation temp. (°C)	37
Normal Low (mg/dl)	2.5
Normal High (mg/dl)	4.5
Linearity Low (mg/dl)	0.17
Linearity High (mg/dl)	23.25
Concentration of standard	See bottle label
Blank with	Reagent blank
Absorbance limit (max.)	2.5
Units	mg/dl

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