

# IMMUNOGLOBULIN G



Cat.No.	Pack name	Packing (Content)
BLT20027	IgG	5 x 25 ml IgG Buffer 1 x 10 ml IgG Antiserum 1 x 1 ml IgG Calibrator



## Intended Use

Quantitative determination of IgG in human serum by turbidimetric immunoassay.

## Diagnostic Implications

IgG is a predominant serum immunoglobulin. The measurement of IgG is important for typing immunodeficiencies and myelomas. Increased levels are found in chronic infections and chronic inflammation. IgG is the only immunoglobulin which crosses the placenta and is therefore of special importance in the infants defense against infection.

## Method

Measurement of antigen-antibody reaction by the end-point method.

## Reagents Provided

### Buffer

Phosphate buffered saline (pH 7.43).  
Polyethylene glycol (5 g/l)  
Sodium azide (0.95 g/l).

### Antiserum

Phosphate buffered saline (pH 7.43).  
Polyclonal goat anti-human IgG (variable).  
Sodium azide (0.95 g/l).

### Calibrator

Pooled human serum, liquid and stabilized. Contains 0.09% sodium azide as preservative.  
Concentration: See bottle label

## Preparation and Stability of Reagents

### Reagent Preparation

Liquid reagents, ready for use.

### Stability and Storage

The reagents are stable until expiry date when kept at 2-8°C. Stability in the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

### Reagents required but not supplied

- Sodium Chloride (9 g/l)
- Calibrators and Controls

12000153

Cat. No.:	Product name	Pack name	Content
BLT20024	IgG CALIBRATOR	IgG CAL H	1 x 1 ml
BLT20034	MULTICONTROL LEVEL 1	MULTICON L1	1 x 1 ml
BLT20035	MULTICONTROL LEVEL 2	MULTICON L1	1 x 1 ml

## Sample collection

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2 - 8°C for 48 hours. If stored for a longer period, the sample should be frozen.

## Automation

Application procedures on clinical chemistry analyzers are available upon request.

## Manual Procedure

Sample/Control: dilute 1:10 in saline 9 g/l  
Reference curve : generate a reference curve by diluting the standard high level 1:10, 1:20, 1:40, 1:80, 1:160 in saline 9 g/l. Use saline 9 g/l as zero point.  
Test : Mix 20 µl diluted standards, controls and samples with 900 µl buffer. Read optical density (OD1) of samples, controls and standards at 340 nm. Add 100 µl of IgG Antiserum. Mix and incubate for 5 minutes at room temperature. Read optical densities (OD2) of samples, controls and standards at 340 nm.

Calculate  $\Delta OD$ 's, plot a standard curve and read the concentration of controls and samples.

## Reference Values

680 - 1445 mg/dl (IFCC)  
This range is given for orientation only. Each laboratory should establish its own reference values.

## Performances

The performance characteristics for the IgG reagents were measured on a clinical chemistry analyzer.

Measuring Range : 0 - 2700 mg/dl

Detection Limit : 90 mg/dl

Hook effect : > 10000 mg/dl

Sensitivity : 0.00014 ABS units/concentration unit

Precision : [%CV]

	Low	Medium	High
Intra-Run	2.50	3.25	4.18
Inter-Run	4.08	1.83	

Accuracy : [mg/dl]

Control	Assigned	Measured
BIO-RAD 1	909(727-1091)	948
BIO-RAD 2	2695(2156-3234)	2732

Specificity :

Monospecific

Interferences : No interference for : Hemoglobin (1000 mg/dl), Na-citrate (1000 mg/dl), Heparin (50 mg/dl), Bilirubin (20 mg/dl), Triglyceride (2500 mg/dl)

Limitations : None

Comparison with nephelometry:  $y = 0.987x - 0.2018$

$$r = 0.9636$$

Stability at 4° C : at least 3 years after production

## Precautions and warnings

- In vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.
- Polyethyleneglycol is not biohazardous
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

QUALITY SYSTEM CERTIFIED  
ISO 9001 ISO 13485




Erba Lachema s.r.o., Karásek 1d, 621 00 Brno, CZ  
e-mail: [diagnostics@erbalachema.com](mailto:diagnostics@erbalachema.com), [www.erbamannheim.com](http://www.erbamannheim.com)  
N/177/16/A/INT Date of Revision: 26.1.2016


## References


Dati, F. et al., Lab. Med. 13, 87 (1989)

### USED SYMBOLS:


 Catalogue No


 CE Mark - Device comply with the Directive 98/79/EC


 Batch Code


 In Vitro Diagnostics

 Expiry Date.

 Consult Instruction for Use

 Manufactured by

 Storage temperature

 Product Name

 Content

