GLUCOSE

Cat. No.	Pack Name	Packaging (Content)
XSYS0012	GLU 440	R1: 10 x 44 ml
XSYS0069	GLU 576 XL-1000	R1: 8 x 72 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Glucose in human serum, plasma and urine.


CLINICAL SIGNIFICANCE

Accurate measurement of glucose in body fluid is important in diagnosis and management of diabetes, hypoglycemia, adrenal dysfunction and various other conditions. High levels of serum glucose may be seen in case of Diabetes mellitus, in patients receiving glucose containing fluids intravenously, during severe stress and in cerebrovascular accidents. Decreased levels of glucose can be due to insulin administration, as a result of insulinoma, inborn errors of carbohydrate metabolism or fasting.

PRINCIPLE

Trinder's method

Glucose in the sample is oxidised to yield gluconic acid and hydrogen peroxide in the presence of Glucose oxidase. The enzyme peroxidase catalyses the oxidative coupling of 4-aminoantipyrine with phenol to yield a coloured quinonemine complex, with absorbance proportional to the concentration of glucose in sample.

REAGENT COMPOSITION

R1

Phosphate buffer	250 mmol/l
Glucose oxidase	> 25 U/ml
Peroxidase	> 2 U/ml
Phenol	5 mmol/l
4-aminoantipyrine	0.5 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

On board stability: min. 30 days if refrigerated (2–10°C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability after addition of a glycolytic inhibitor (Fluoride, monoiodoacetate, mannose): ⁴

2 days	at 20–25°C
7 days	at 4–8°C

Stability in serum (separated from cellular contents, hemolysis free) without adding a glycolytic inhibitor: ^{2,5}

8 hours	at 25°C
72 hours	at 4°C

Stability of glucose in urine:

24 hours at 4–8°C

For the determination in urine dilute the sample using redistilled water in 1 + 10 ratio. Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to ID-MS.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

mg/dl x 0.056 = mmol/l

EXPECTED VALUES ²

Serum:

Glucose Fasting:

Cord:	45 – 96 mg/dl
Newborn, 1 d:	40 – 60 mg/dl
Newborn, >1 d:	50 – 80 mg/dl
Child:	60 – 100 mg/dl
Adult:	74 – 100 mg/dl
>60 y:	82 – 115 mg/dl
>90 y:	75 – 121 mg/dl

Glucose 2 h Postprandial: <120 mg/dl

WB (Hep) Adult: 65 – 95 mg/dl

Urine: 1 – 15 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	2.34 mg/dl
Linearity:	450 mg/dl
Measuring range:	2.34 – 450 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	153.75	1.61	1.05
Sample 2	239.64	3.93	1.66

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	97.14	1.61	1.64
Sample 2	258.39	3.04	1.19

COMPARISON

A comparison between XL-Systems Glucose (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.000 x - 0.714 mg/dl

r = 0.999

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 7.5 g/l, bilirubin up to 30 mg/dl, triglycerides up to 750 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent R1 contains 0.74% NaOH and is classified as irritant and also contains less than 0.1% sodium azide which is classified as very toxic and dangerous substance for the environment.

Warning

Hazard statements (H):

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statements (P):

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	GLU	GLU	GLU	GLU	GLU	GLU
Test Code	23	23	23	23	23	23
Report Name	Glucose	Glucose	Glucose	Glucose	Glucose	Glucose
Unit	mg/dl	mg/dl	mg/dl	mg/dl	mg/dl	mg/dl
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	505	505	505	505	505	505
Wavelength-Secondary	660	660	660	660	660	660
Assay type	1-Point	1-Point	1-Point	1-Point	1-Point	1-Point
Curve type	Linear	Linear	Linear	Linear	Linear	Linear
M1 Start	0	0	0	0	0	0
M1 End	0	0	0	0	0	0
M2 Start	32	34	48	60	29	32
M2 End	34	36	51	63	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	2.2	2.2	2.2	2.2	2.2	2.2
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	2.34	2.34	2.34	2.34	2.34	2.34
Technical Maximum	450	450	450	450	450	450
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	0	0	0	0	0	0
Reagent Abs Max	0.2	0.2	0.2	0.2	0.2	0.2
Auto Rerun	No	No	No	No	No	No
Total Reagents	1	1	1	1	1	1
Reagent R1	GLU R1	GLU R1	GLU R1	GLU R1	GLU R1	GLU R1
Reagent R2	NA	NA	NA	NA	NA	NA
Reagent R3	NA	NA	NA	NA	NA	NA

Test Volumes						
Test	GLU	GLU	GLU	GLU	GLU	GLU
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes						
Normal	2	2	2	2	2	2
Dilution Ratio	1	1	1	1	1	1
Increase	4	4	4	4	4	4
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	2	2
Dilution Ratio	5	5	5	5	5	5
Standard volume	2	2	2	2	2	2
Reagent Volumes and Stirrer speed						
RGT-1 Volume	200	200	200	200	200	200
R1 Stirrer Speed	High	High	NA	High	High	High
RGT-2 Volume	0	0	0	0	0	0
R2 Stirrer Speed	NA	NA	NA	NA	NA	NA
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA

Reference Ranges						
Test	GLU	GLU	GLU	GLU	GLU	GLU
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	74	74	74	74	74	74
Normal-Upper Limit	100	100	100	100	100	100
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	74	74	74	74	74	74
Normal-Upper Limit	100	100	100	100	100	100
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA

Revision Number						
Revision	<A-100- GLU-1 20.08.2013>	<A-200- GLU-1 20.08.2013>	<A-300/600- GLU-1 20.08.2013>	<A-640- GLU-1 20.08.2013>	<A-1000- GLU-1 20.08.2013>	<A-180- GLU-1 12.12.2013>

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1. Thomas L.: Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 131 - 7.
2. Tietz N. W., (Ed.), Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
3. Barham, D., Trinder, P.: An improved color reagent for the determination of blood glucose by the oxidase system. Analyst, 1972, 97; 142 - 5.
4. Guder WG, Zawta B et al. The quality of Diagnostic Samples. 1st ed. Darmstadt: GIT verlag; 2001; p.30-1.
5. Snacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chemi 2002; 48:436-72.

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