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Literatúra

- Rifai, N., Bachorik, PS., Albers, JJ.: Lipids, lipoproteins and apolipoproteins. In: Burtis, CA., Ashwood, ER., editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company, 1999, p. 809 - 61.
- Cole, TG., Klotzsch, SG., McNamara, J.: Measurement of triglyceride concentration. In: Rifai, N., Warnick, GR., Dominiczak, MH., eds. Handbook of lipoprotein testing. Washington: AACC Press, 1997, p. 115 - 503.
- Recommendation of the Societies on Coronary Prevention. Prevention of coronary heart disease in clinical practise. Eur Heart J., 1998, 19: p. 1434 - 503.

Dátum poslednej revízie: 15. 4. 2011

IVD


TRIACYLGLYCEROLS LIQUID 250 S (TG L 250 S)

Cat. No. 10003267

Store at (+2 to +8) °C

Reagent for enzymatic photometric determination of triacylglycerols in serum or plasma.

Principle

Reagents

R1 Reagent

Good's buffer (pH 7.2)	1 x 250 ml
4-chlorophenol	50 mmol/l
Mg ²⁺	4 mmol/l
ATP	15 mmol/l
Glycerolkinase	2 mmol/l
Peroxidase	≥ 6.7 µkat/l
Lipoproteinlipase	≥ 33 µkat/l
Glycerol-3-phosphate-oxidase	≥ 33 µkat/l
4-aminoantipyrine	≥ 8.3 µkat/l
	0.5 mmol/l

R2 Standard

Triacylglycerols	1 x 3 ml
	2.3 mmol/l

Reaction mixture

Good's buffer (pH 7.2)	49.5 mmol/l
4-chlorophenol	3.96 mmol/l
Mg ²⁺	14.85 mmol/l
ATP	1.98 mmol/l
Glycerolkinase	≥ 6.6 µkat/l
Peroxidase	≥ 32.7 µkat/l
Lipoproteinlipase	≥ 32.7 µkat/l
Glycerol-3-phosphate-oxidase	≥ 8.2 µkat/l
4-aminoantipyrine	0.495 mmol/l

Preparation and stability of working reagents

Reagents R1 and R2 are liquid, ready to use.

If stored at (+2 to +8) °C, reagent is stable until expiry date, that is stated on the package.

After opening, reagent is stable until expiry date at (+2 to +8) °C if stored at appropriate conditions, in the dark place, closed carefully and without any contamination.

Samples

Non – haemolytic serum, plasma (heparin, EDTA)

Stability of triacylglycerols in serum and plasma is 10 days at (+4 to +8) °C.

Calibration

For the calibration, it is recommended to use Lyonorm Calibrator, Cat. No. 10003200.

Quality control

For quality control, it is recommended to use:

Lyonorm HUM N, Cat. No. 10003204

Lyonorm HUM P, Cat. No. 10003206

Procedure

Wavelength	500 nm (546 nm)
Cuvette	1 cm
Temperature	(+20 to +25) °C, +37 °C
Serum/reaction mixture ratio	1/101

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

	Reagent blank	Standard (calibrator)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	–	–	0.01 ml
Standard (calibrator)	–	0.01 ml	–
Distilled water	0.01 ml	–	–

Mix and incubate 10 min. at 37 °C or 20 min. at 20–25 °C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank. The coloration is stable during one hour.

Applications for automatic analysers will be supplied on request.

Calculation

$$\text{Triacylglycerols (mmol/l)} = \frac{\Delta A_{sam}}{\Delta A_{st}} \times C_{st}$$

C_{st} = standard (calibrator) concentration

Unit conversion

1 mmol/l = 0.0113 mg/dl

Reference values

fS Triacylglycerols (mmol/l)

up to 1.92

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

Performance data

Linearity: up to 11.4 mmol/l

Limit of detection: 0.03 mmol/l

Low limit of quantification: 0.1 mmol/l

Working range: 0.1 – 11.4 mmol/l

Precision (at 37 °C)

INTRA-ASSAY n = 20	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	1.13	0.032	2.28
Sample 2	1.42	0.020	1.48
Sample 3	2.33	0.027	1.15
INTER-ASSAY n = 10	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	1.16	0.036	3.12
Sample 2	1.42	0.033	2.31
Sample 3	2.33	0.032	1.36

Accuracy

Bias

–3.60% (1.39 mmol/l)

–2.94% (2.38 mmol/l)

Correlation: with commercial available method. Linear regression

n = 97, r = 0.995, y = 0.9397 x + 0.0312 mmol/l

Interference

Following substances do not interfere:

Ascorbic acid up to 0.34 mmol/l, bilirubin up to 684 µmol/l, haemoglobin up to 2.5 g/l.

Notes

- The measurement is not influenced by occasionally occurring color changes, as long as the absorbance of the reagent is < 0.3 (at 546 nm).
- If the triacylglycerols concentration exceeds 11.4 mmol/l, dilute the sample by saline (0.9% NaCl) in ratio 1 + 4 (result x 5).

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

References

- Rifai, N., Bachorik, PS., Albers, JJ.: Lipids, lipoproteins and apolipoproteins. In: Burtis, CA., Ashwood, ER., editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company, 1999, p. 809 - 61.
- Cole, TG., Klotzsch, SG., McNamara, J.: Measurement of triglyceride concentration. In: Rifai, N., Warnick, GR., Dominiczak, MH., eds. Handbook of lipoprotein testing. Washington: AACC Press, 1997, p. 115 - 503.
- Recommendation of the Societes on Coronary Prevention. Prevention of coronary heart disease in clinical practise. Eur Heart J., 1998, 19: p. 1434 - 503.

Date of last revision: 15. 4. 2011