**APOLIPOPROTEIN A1
READY TO USE****(IT-Apo A1 DIL 1x15)**

Kat. č. 10006879

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Literatúra

- Marcovina, S. M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B. Clin. Chem. 37/10, (1991) 1676-1682
- Marcovina, S. M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B.III. Comparability of apolipoproteins A1 values by use of international reference material Clin. Chem. 39, (1993) 773-781
- Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Zima, T.: Laboratórna diagnostika, Galén 2002

Dátum poslednej revízie: 7. 9. 2010

**APOLIPOPROTEIN A1
READY TO USE****(IT-Apo A1 DIL 1x15)**

Cat. No. 10006879

Store at (+2 to +8) °C

Reagent kit for immunoturbidimetric determination of apolipoprotein A1 in human serum and plasma.

Principle

The human serum (plasma) sample reacts upon a specific antibody for human apolipoprotein A1 and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the apolipoprotein A1 sample concentration.

Reagents

R1 Buffer	(1x125 ml)
TRIS pH 8	200 mmol/l
PEG	56 g/l
Sodium chloride	50 mmol/l
Sodium azide	0.9 g/l
R2 Antiserum	(1x15 ml)
Sterile delipidated goat serum against human Apo A1	titre +/- 0.5 g/l
HEPES buffer pH 7.4	50 mmol/l
EDTA	9 mmol/l
Sodium azide	0.9 g/l

Reaction mixture

TRIS pH 8	165.66 mmol/l
HEPES buffer pH 7.4	8.28 mmol/l
PEG	46.35 g/l
Sodium chloride	41.4 mmol/l
EDTA	1.49 mmol/l
Sodium azide	0.89 g/l

Preparation and stability of working reagents

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25) °C and reagent R2 is stable until expiry date at (+2 to +8) °C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Samples

Serum or plasma (heparin, NaF, EDTA, sodium citrate), samples should be fresh and not hemolysed.

Prior the analysis, the serum can be stored a week at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

Calibration

For calibration use APO standard (Cat. No. 10003240).

For photometers working with a 5 points calibration curves, prepare a range of 5 reference samples by successive 1:1 dilutions of the APO standard by saline. If, for example, the concentration of apolipoprotein A1 in APO standard is 3.0 g/l, prepare the following concentrations: 0.19-0.38-0.75-1.5-3.0 g/l (dilution in ratio 1:16, 1:8, 1:4, 1:2, undiluted standard).

Plot the calibration curve.

Dilute samples with apolipoprotein A1 concentration higher than the highest standard by saline 10x and repeat analysis. The results multiply by dilution.

Quality control

For quality control use Lyonorm LIPID HUM N (Cat.No.10003286), Lyonorm LIPID HUM P (Cat.No. 10003287), Bio-Rad Liquicheck™ Immunology Control or other quality control sera with apolipoprotein A1 concentration values determined by immunoturbidimetric method.

Procedure

Wavelength	600 nm
Cuvette	1 cm
Temperature	37 °C
Serum/reaction mixture ratio	1/151

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

	Sample	Standard	Reagent blank
Standard	—	0.002 ml	—
Sample	0.002 ml	—	—
Distilled water	—	—	0.002 ml
Reagent 1	0.250 ml	0.250 ml	0.250 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A1.			
Reagent 2	0.050 ml	0.050 ml	0.050 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A2.			

Applications for automatic analysers can be supplied for request.

Calculation

Calculate the difference A2-A1 for standards and samples.

$$\text{ApoA1(g/l)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagent blank}}}{(A_2 - A_1)_{\text{standard}} - (A_2 - A_1)_{\text{reagent blank}}} \times C_{\text{st}}$$

C_{st} ...concentration of standard

Reference values

fS Apo A1 (g/l)

male 0,94 – 1,78

female 1,01 – 1,99

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

Performance characteristics**Limit of detection:** 0.01 g/l**Hook effect:** from 5.0 g/l

Specificity: the antiserum shows the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

Precision (at 37 °C)

INTRA-ASSAY	Mean (g/l)	SD (g/l)	CV (%)
n = 10			
Sample 1	0.76	0.012	1.6
Sample 2	1.50	0.031	2.1
INTER-ASSAY	Mean (g/l)	SD (g/l)	CV (%)
n = 10			
Sample 1	0.76	0.019	2.5
Sample 2	1.52	0.038	2.5

Comparison: with commercial available test (immunoturbidimetry). Results: n = 29, r = 0.977, y = 1.11x – 0.147 g/l

Interferences: Triglycerides to < 22.6 mmol/l, haemoglobin to < 200 µmol/l, bilirubin to < 500 µmol/l do not interfere with the measurements (< 10% or 2 SD).

Calibration stability: on the most often automatic analysers it is > 4 weeks.

Health protection

For in vitro diagnostic use. To be handled by entitled and professionals educated person.

Reagent of the kit are not classified like dangerous.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

Products containing sodium azide have to be handled with care, avoid ingestion and contact with skin and mucous membrane.

First aid

In case of an accidental ingestion wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

References

1. Marcovina, S. M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B. Clin. Chem. 37/10, (1991) 1676-1682
2. Marcovina, S. M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B.III. Comparability of apolipoproteins A1 values by use of international reference material Clin. Chem. 39, (1993) 773-781
3. Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
4. Zima, T.: Laboratorní diagnostika, Galén 2002

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- Marcovina, S.M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B. Clin. Chem. 37/10, (1991) 1676-1682
- Marcovina, S.M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B.III. Comparability of apolipoproteins A1 values by use of international reference material Clin. Chem. 39, (1993) 773-781
- Tietz, N.W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Zima, T.: Laboratorní diagnostika, Galén 2002

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