**REVMATOIDNÍ FAKTOR
READY TO USE****(IC-RF DIL 2x10)**

Kat. č. 10003301

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**RHEUMATOID FACTOR
READY TO USE****(IC-RF DIL 2x10)**

Cat. No. 10003301

Store at (+2 to +8) °C

Reagent kit for immunocolorimetric determination of rheumatoid factor in human serum and plasma.

Principle

The human serum (plasma) sample reacts upon colloidal gold particles coated with human IgG. In the presence of IgM rheumatoid factor, the particles agglutinate, which induces a red shift in a visible spectrum of a colloid. This induces an increase in optical density at 600 nm, which is directly proportional to the rheumatoid factor concentration in the sample.

Reagents

R1 Pufr	(1 x 40 ml)
TRIS buffer pH 8	50 mmol/l
PEG	68 g/l
Sodium chloride	50 mmol/l
EDTA	10 mmol/l
Sodium azide	< 1 g/l
R2 RF reagent	(2 x 10 ml)
Human IgG coated on colloidal gold particles	titre +/- 60 mg/l
Veronal buffer pH 8	50 mmol/l
Sodium chloride	50 mmol/l
EDTA	9 mmol/l
Sodium azide	< 1 g/l

Reaction mixture

TRIS pufr pH 8	31 mmol/l
Veronal buffer pH 8	18.4 mmol/l
PEG	42 g/l
Sodium chloride	49.4 mmol/l
EDTA	9.5 mmol/l
Sodium azide	< 1 g/l

Preparation and stability of working reagents

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25) °C and reagent R2 is stable until expiry date at (+2 to +8) °C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Samples

Serum or plasma (heparin, NaF, EDTA, sodium citrate), samples should be fresh and not hemolysed.

Prior the analysis, the serum can be stored a week at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

Calibration

For calibration use RF standard (Cat. No. 10003245).

For photometers working with a 5 points calibration curves, prepare a range of 5 reference samples by successive 1:1 dilutions of the RF standard by saline. If, for example, the concentration of rheumatoid factor in RF standard is 140 IU/ml, prepare the following concentrations: 8.8-17.5-35-70-140 IU/ml (dilution in ratio 1:16, 1:8, 1:4, 1:2, undiluted standard).

Plot the calibration curve.

Dilute samples with rheumatoid factor concentration higher than the highest standard by saline 10x and repeat analysis. The results multiply by dilution.

Quality control

For quality control use Bio-Rad Liquichek™ Rheumatoid Factor Control or other quality control sera with rheumatoid factor concentration values determined by immunocolorimetric method.

Procedure

Wavelength	600 nm
Cuvette	1 cm
Temperature	37 °C
Serum/reaction mixture ratio	1/86

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

	Sample	Standard	Reagent blank
Standard	–	0.003 ml	–
Sample	0.003 ml	–	–
Distilled water	–	–	0.003 ml
Reagent 1	0.160 ml	0.160 ml	0.160 ml
Reagent 2	0.095 ml	0.095 ml	0.095 ml
Mix R1 and R2, add sample and read the absorbance A1 within 30 seconds. After 5-10 minutes incubation at 37 °C read the absorbance A2.			

Applications for automatic analysers can be supplied for request.

Calculation

Calculate the difference A2–A1 for standards and samples.

$$RF \text{ (IU/ml)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagent blank}}}{(A_2 - A_1)_{\text{st}} - (A_2 - A_1)_{\text{reagent blank}}} \times C_{\text{st}}$$

C_{st} ...concentration of standard

Reference values

fS rheumatoid factor (IU/l)

< 20

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

Performance characteristics

Limit of detection: 2.1 IU/l

Hook effect: from 1000 IU/l

Specificity: the antiserum shows the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

Precision (at 37 °C)

INTRA-ASSAY n = 10	Mean (IU/l)	SD (IU/l)	CV (%)
Sample 1	14.0	0.40	3.2
Sample 2	24.0	0.80	3.4
INTER-ASSAY n = 10	Mean (IU/l)	SD (IU/l)	CV (%)
Sample 1	13.0	0.56	4.3
Sample 2	19.6	0.69	3.5

Comparison: with commercial available test (latex immunoturbidimetry).

Results: n = 55, r = 0.985, y = 0.94x + 2.51 IU/ml

Interferences: Triglycerides to < 4.5 mmol/l, hemoglobin to < 200 µmol/l, bilirubin to < 500 µmol/l do not interfere with the measurements (< 10 % or 2 SD).

Calibration stability: on the most often automatic analysers it is > 4 weeks.

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

Products containing sodium azide have to be handled with care, avoid ingestion and contact with skin and mucous membrane.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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