IVD

Činidlo pro test PYR

CZ

Kat. č.: MLT00023

Pro mikrobiologii

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Cat. No.: MLT00023

For microbiology

Reagent for PYR test is the colour forming reagent for performing of pyrrolidonylarylamidase test (PYRAtest strips and/or test PYR in MIKROLATEST® kits). The working instructions for employing the reagent are a part of the working instructions of the respective MIKROLATEST® kit.

Reagent for PYR test contains:

Reagent for PYR test in dropping bottle 18 ml
Information leaflet 1 pc

Storage, expiration:

Reagent for PYR test should be stored at the temperature of (+2 až +8) °C until the expiry date indicated on the package. **Note:** in the solution of Reagent for PYR may after prolonged storage to form a suspension which is not harmful.

Warning:

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent for PYR contain harmful 2-methoxyethanol (< 18 %), which is toxic for fertility (category 1B) and Triton X-100 (< 3.5 %).

Danger

Hazard statement:

H318 Causes serious eye damage.

H360FD May damage fertility. May damage the unborn child.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P308 + P313 If exposed or concerned: Get medical advice/lattention.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If swallowed it is necessary to drink about 0.5 l of lukewarm water. On skin or eye contact rinse the place with jet of tap water.

If necessary, consult a physician.

Waste disposal:

Contingent rest of the reagent should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

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