CRP CALIBRATOR

Lot: XXXXXXX

Cat. No.	Pack Name	Packaging (Content)
XSYS0053	CRP CAL SH	1 x 1 ml

INTENDED USE

Preparation of reference curves for quantitative immunochemical determination of CRP.

COMPOSITION

A dilution of human plasma and pleural fluid containing high levels of CRP with phosphate buffered saline.

Liquid stabilised.

Ready to use.

STORAGE AND STABILITY

The expiry date of the product at 2–8 °C is listed on the label.

After first opening the container, the serum can be used for 6 weeks if stored tightly closed at 2–8 °C after use.

Do not freeze.

PRECAUTIONS AND WARNING

- For *in Vitro* Diagnostic use.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.

ASSIGNED VALUES

The value for CRP has been assigned to ERM-DA474/IFCC

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