# HbA1c CONTROL

Lot: XXXXXXX

Cat. No.	Pack Name	Packaging (Content)
XSYS0055	HBA1C CON L	1 x 0.25 ml

## INTENDED USE

HbA1c liquid control is a human-based control. The HbA1c concentration in HbA1c Control Low is normal.

## STORAGE AND STABILITY

The controls both unopened and opened must be stored at 2–8 °C, protected from light and heat.

Unopened and opened: until the end of the indicated month of expiry if contamination and evaporation are avoided after having opened the vials.

Proper storage and handling of this product must be observed.

## WARNINGS AND PRECAUTIONS

Reagent of the kit is not classified like dangerous but contains less than 0.1 % sodium azide classified as very toxic and dangerous substance for the environment.

Each individual blood donation used for production of HbA1c control was found to be non-reactive when tested with approved methods for HBsAg, anti-HIV-1+2 and anti-HCV. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the control with the same precautions used for patient specimens.

## PREPARATION

HbA1c control is ready to use. Controls must be treated the same way as patient samples. Please refer to the package insert of the reagent HbA1c, Cat. No. XSYS0054

## CONTROL PREPARATION:

Hemolyzing Solution (R3) 500 µl  
Control 10 µl  
Mix and allow to stand for 5 minutes or until complete lysis is apparent.

## ASSIGNED VALUES

Target values according to IFCC (mmol/mol):

	Lot No.	Expiry Date		Targer value (mmol/mol)	Range (mmol/mol)	1 SD (mmol/mol)
HbA1c Control Low	xxxxxxx	mm/yyyy	3-component system	xx,x	xx,x – xx,x	x,xx
			2-component system	xx,x	xx,x – xx,x	x,xx

Target values according to DCCT/NGSP (%):

	Lot No.	Expiry Date		Targer value (%)	Range (%)	1 SD (%)
HbA1c Control Low	xxxxxxx	mm/yyyy	3-component system	x,xx	x,xx – x,xx	x,xx
			2-component system	x,xx	x,xx – x,xx	x,xx

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