Nr kat.: MLT00080

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The kit is designed to test antimicrobial susceptibility of staphylococci on the basis of bacterial growth at breakpoint concentrations. The kit contains 60 detection strips. Breakpoint concentrations and interpretation of sensitivities were determined on the basis of the European standard EUCAST ([www.eucast.org](http://www.eucast.org)) dated from 5.1.2011. The test is based on the rehydration of antibiotics in the wells with Suspension medium for SENSILATEST and addition of bacterial suspension. The results are read visually after 18-24 hours of incubation.

**The kit contains:**

- 10 plates for examination
- A lid
- 10 pc of PE bags
- Record sheet
- Storage bag
- 1 pc of Suspension medium for SENSILATEST in dehydrated form – Attention: It is necessary to prepare it before testing!

To avoid preparation of suspension media from dehydrated form, it is possible to order ready to use media. Ready to use media are available as item MLT00082 - Suspension medium for SENSILATEST. Each kit contains 20 glass tubes. One tube contains 4.4 ml of suspension media, which allows one examination on any type of SENSILATEST strip.

**Storage and expiration of the kit:**

It is recommended to store the kit at (+15 to +25) °C. The date of expiration is indicated on each package. After the aluminium package is opened, it is possible to insert any unused strips into a storage bag with desiccant at the room temperature for a period of maximum 8 days. Mark the storage bag with a type of stored strips (e.g. G-I, G-II etc.) to avoid mistake in their further usage. **Don't leave opened strips unprotected!!! Exposition to open air leads to antibiotic activity failure!!!** In case that laboratory is not able to keep kits at recommended room temperature (not exceeding 25 °C), it is better to keep kits in refrigerator (+2 to +8) °C. Then it is necessary to leave plate at room temperature for 2 hours before the use to avoid water condensation.

**Storage and expiration of Suspension medium for SENSILATEST:**

It can be stored in powder form at (+15 to +25) °C till the date indicated on the package. The expiration of the Suspension medium for SENSILATEST is 6 months at (+15 to +25) °C after it was prepared into a liquid form. Do not store the prepared medium at lower temperatures as it can lead to formation of gel in the media. If such a situation occurs, adjust the medium to room temperature and shake it thoroughly before the use.

**Material required to perform a test, not included in the kit:**

- Sterile physiological solution
- Sterile tubes
- An stepper pipette for dosage of 100 µl
- A pipette for dosage of 40 µl
- Densitometer (e.g. Densi-La-Meter II, Erba Lachema)
- Incubator 37 °C
- Regular microbiological laboratory equipment (loops, marker, burner, etc.)

**Material required to prepare Suspension medium for SENSILATEST, not included in the kit:**

- A mixing stick or a stirrer
- Autoclave

**Caution:**

**The kit is for professional use only!**  
**Respect the rules for work with infectious material!**

## Instructions for Use

**Caution:**

**Prepare the Suspension medium for SENSILATEST before you start testing!**

To avoid preparation of suspension media from dehydrated form, it is possible to order ready to use media. Ready to use media are available as item MLT00082 - Suspension medium for SENSILATEST. Each kit contains 20 glass tubes. One tube contains 4.4 ml of suspension media, which allows one examination on any type of SENSILATEST strip.

**Preparation of Suspension medium for SENSILATEST:**

One pack is sufficient to prepare 0.5 L of Suspension medium for SENSILATEST. Add 0.5 L of preheated distilled water and stir it. Autoclave the prepared suspension for 30 minutes at 121 °C. Cool it to room temperature and stir it again under sterile conditions. A small precipitate at the bottom does not influence the function of the medium. Check the pH of the medium before the use. pH of the medium must be 7.5 (± 0.1). Adjust pH if necessary.

**Instructions to use the kit:**

- 1) Prepare a tube with 2 - 5 ml of physiological solution of pH 5.5 – 6.5.
- 2) Prepare a tube with Suspension medium for SENSILATEST optimally preheated to 37 °C.
- 3) Remove few colonies from 18 – 24 hour culture on blood agar and prepare a bacterial suspension of density of 0.5 on McF scale in physiological solution.  
**Overestimation of suspension density is smaller error than its underestimation.**
- 4) Inoculate **40 µl** of bacterial suspension into 4.4 ml of Suspension medium for SENSILATEST and homogenise well.

**Inoculation:**

Remove a plate from aluminium bag. Remove required number of strips for examination. Do not remove the strips from the package earlier than 30 minutes before you start the work. Insert any unused strips into a ZIP log bag. Transfer a desiccant bag from the original package into a ZIP log bag and close it for later use. Mark the storage bag with a type of stored strips (e.g. G-I, G-II etc.) to avoid mistake in their further usage. **Use the strips stored in this way within 8 days! Don't leave opened strips unprotected!!! Exposition to open air leads to antibiotic activity failure!!!**

Record numbers of the examined strains on the corresponding strips. Mark the frame with a type of strip (e.g. G-I, G-II etc.) to avoid mistake in reading results after incubation.

Inoculate **100 µl** of suspension in Suspension medium for SENSILATEST into each well of the strip with the help of automatic stepper pipette.

**Incubation:**

Insert the inoculated plate into a PE bag. Fold the open end of the bag under the plate to prevent evaporation during the incubation. Incubate the plate at 37 °C for 18 – 24 hours.

**Evaluation:**

Remove the plate from the PE bag. To read the growth in the microwells, choose a way which is the most convenient for you:

- 1) Use the lid with imprinted antibiotics sequence.
- 2) Read against a grey background.
- 3) Read against natural or artificial dispersed light – you can detect even limited growth in a form of grains or light turbidity.

**Please read with attention:**

You must see a growth in a control well! If the growth is not present, the test **MUST NOT** be evaluated! However, the growth control is not used to compare the presence of the growth in individual wells with antibiotics. Every change from a total transparency is evaluated as growth, including small grains of growth in a well, light turbidity or partial growth covering only a part of a well. Beware to differentiate grains of growth from media bubbles. Record the results.

| H              | G                             | F            | E             | D          | C            | B          | A          |
|----------------|-------------------------------|--------------|---------------|------------|--------------|------------|------------|
| Growth control | Trimetoprim / sulfamethoxazol |              | Tigecycline   | Linezolid  | Erythromycin |            | Vancomycin |
| GC             | T/S<br>2/38                   | T/S<br>4/76  | TGC<br>0,5    | LIZ<br>4   | ERY<br>1     | ERY<br>2   | VAN<br>2   |
| Ciprofloxacin  | Cefoxitin                     | Fusidic acid | Rifampicin    |            | Clindamicin  |            | Gentamicin |
| CIP<br>1       | COX<br>4                      | FUS<br>1     | RAM<br>0,0625 | RAM<br>0,5 | CLI<br>0,25  | CLI<br>0,5 | GEN<br>1   |

**Interpretation table:**

|                               |     |   |   |   |   |   |   |   |   |   | Results: +/- | Evaluation (S, I, R) |
|-------------------------------|-----|---|---|---|---|---|---|---|---|---|--------------|----------------------|
| Growth control                | GC  | + |   |   |   |   |   |   |   |   |              |                      |
| Cefoxitin                     | COX | - |   | S |   |   |   | + |   | R | x            |                      |
| Clindamicin                   | CLI | - | - | S | + | - | I | + | + | R |              |                      |
| Ciprofloxacin                 | CIP | - |   | S |   |   |   | + |   | R | x            |                      |
| Erythromycin                  | ERY | - | - | S | + | - | I | + | + | R |              |                      |
| Gentamicin                    | GEN | - |   | S |   |   |   | + |   | R | x            |                      |
| Fusidic acid                  | FUS | - |   | S |   |   |   | + |   | R | x            |                      |
| Linezolid                     | LIZ | - |   | S |   |   |   | + |   | R | x            |                      |
| Rifampicin                    | RAM | - | - | S | + | - | I | + | + | R |              |                      |
| Tigecycline                   | TGC | - |   | S |   |   |   | + |   | R | x            |                      |
| Trimetoprim / sulfamethoxazol | T/S | - | - | S | + | - | I | + | + | R |              |                      |
| Vancomycin                    | VAN | - |   | S |   |   |   | + |   | R | x            |                      |

If there is a growth at a higher concentration of an antibiotic and no growth at a lower concentrations of the same antibiotic, do not evaluate susceptibility to this antibiotic. This combination shows an error in the procedure!

In case of **ERY** resistant along with **CLI** sensitive or intermediate, it is necessary to determinate “induction“ of **MLSb** (macrolides-lincosamides-streptogramins B) resistance with disks ERY and CLI and correct the result for CLI.

Resistance to **COX** in **S. aureus** indicates possibility of MRSA. MRSA is necessary to be confirmed by other methods – growth on chromogenic agar, latex agglutination (PBP2a) or PCR.

Resistance to **VAN** is necessary to be confirmed by other methods (gradient MIC, e.g. E-test, or PCR).

Depending on national or laboratory standards, other criteria described in EUCAST Expert Rules in Antimicrobial Susceptibility Testing can be used.

**Quality control:**

Control strains can be used for internal testing of functionality of the antibiotics in the laboratory. Please refer to the acceptable limits for quality control strains according to EUCAST.

| CCM 4223 <i>Staphylococcus aureus</i> (ATCC 29213) |               |               |                 |              |               |              |            |
|----------------------------------------------------|---------------|---------------|-----------------|--------------|---------------|--------------|------------|
| GC+                                                | T/S –<br>2/38 | T/S –<br>4/76 | TGC –<br>0,5    | LIZ –<br>4   | ERY –<br>1    | ERY –<br>2   | VAN –<br>2 |
| CIP –<br>1                                         | COX –<br>4    | FUS –<br>1    | RAM –<br>0,0625 | RAM –<br>0,5 | CLI –<br>0,25 | CLI –<br>0,5 | GEN –<br>1 |

+ growth – no growth

**Other recommended control strains:**

CCM 3954 *Escherichia coli* (ATCC 25922)

CCM 3955 *Pseudomonas aeruginosa* (ATCC 27853)

ATCC – American Type Culture Collection

CCM – Czech Collection of Microorganisms

**Health protection:**

Components of the kit do not contain any dangerous substances.

**Disposal of the used material:**

Insert the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration.

Put paper packaging waste to recycling.

Date of revision: 0. 6. 2019

**USED SYMBOLS**

Catalogue number

In vitro diagnostics

Manufacturer

See instruction for use

Lot number

Storage temperature

Expiry date